



1. Installing the Secure Holter Transfer Program:

- Visit www.HolterLabs.com and click on 'Install Holter Software' link at the top of the webpage.
- A new window will open showing the 'Secure Transfer Application' install icon.
- Click on "Launch Now" to install the program onto your computer.

2. Sending us the Patient Holter Data:

- Remove the 'Memory Card' from the Holter Recorder after the patient has worn the device.
- Place the 'Memory Card' into the **Card Reader** on your computer.
(The Card Reader was supplied our office, just plug the card reader into any USB slot.)
- Run the 'Holter Application' by double clicking the icon that was installed on your desktop.
- Once opened, logon using the credentials given to you by our staff and click 'Login'.
Username: Last name, first initial (ie: Alex Smith = smitha)
Password: password101 (you may change the password anytime; see below)
* Username and password are all in lower case.
- Enter all of the information requested by the program. (ie: Physician, Patients Name, DOB, etc)
- Once on 'Step 4' click on "Upload Files" to send us the patients' data from the memory card.
- In about 5-10 mins; a pop-up message will appear stating that the data was transferred successfully.
- FAX** the 'Holter Order Form' (located in the info binder) to: **888-821-4677**.



3. Viewing the completed Holter reports:

- Logon to the program with the credentials given to you by our staff if you have not done so already.
- Double-click on the patients name in the "Reports" window, and wait for the report to open.
- Screen shot shown below of the Holter Application.

Changing your password:

You may change your default password by visiting www.HolterLabs.com and clicking on the 'Doctor Sign In' link.

Logon with your current credentials and enter in a new password in the "change password" section toward the bottom of the page.

Software & Data Transfer
HIPPA Compliant



The screenshot shows the 'New Upload' window of the Holter Labs Cardiac Monitoring software. The window has a dark grey background with the Holter Labs logo at the top left. On the right side, it says 'Hello Dr. [Name]' and 'Logout'. The main area is divided into two sections: 'Reports' on the left and 'New Upload' on the right. The 'Reports' section lists several patient names and dates, such as 'Acosta, Kathryn - 3-24-08.pdf'. The 'New Upload' section has a progress bar with four steps, where step 2 is highlighted. Below the progress bar, there are several input fields for 'Patient Information': 'Patient Name:', 'Patient DOB:', 'Date of Service:', 'Time of Service:', and 'Comments:'. Each of these fields has a red 'Required' label next to it. At the bottom of the 'New Upload' section, there are two buttons: 'Previous Step' and 'Next Step'. At the very bottom of the window, there is a red bar with the text 'SECURE TRANSFER APPLICATION' and 'Problems? Contact us at 1.888.821.4667 or support@holterlabs.com'.